Rapid response to high cut-off haemodialysis and bortezomib therapy in a patient with acute renal failure and plasma cells dyscrasia

Nefrologia 2014;34(5):688

doi:10.3265/Nefrologia.pre2014.Jun.12108

To the Editor,

Renal involvement is a complication of plasma cell dyscrasias. Electrochemical properties of abnormal light chains responsible for histological renal pictures. Bortezomib is the first proteasome inhibitor approved for treatment of multiple myeloma and amyloidosis. It prevents the activation of NF-kB that controls the genes encoding IL-6, TNF-alpha and other cytokines and growth factors. Our experience: a 29-yearold woman, suffering from ulcerative colitis successfully treated sulfasalazine, without renal disease history, was admitted to hospital for acute kidney injury (creatinine 10mg/dL) and anaemia without signs of thrombotic microangiopathy. Serum protein electrophoresis showed a beta-2-monoclonal peak; the serum immunofixation was positive for λ light chains; k- and λ -free light chains (FLC) levels were 37mg/L and 1750mg/L respectively with k/λ ratio=0.02. Microbiological, coagulation other immunological investigations unremarkable. The patient started haemodialysis treatment. Renal biopsy revealed a cast-nephropathy picture, negative Congo-red staining, without glomerular deposits; the immunofluorescence showed k light chains but not λ light chains in tubular basement membranes. Abdominal fat biopsy was positive for AL amyloid. Bone marrow biopsy showed 10% of plasma cells infiltration without morphological or cytometric clonality markers; total body CT scan was negative for bone lesions. Because of plasma cell dyscrasia and severe renal

involvement, even in absence of a diagnostic definition, was performed a combination treatment of direct removal of FLCs and chemotherapy. Following Hutchison's studies,1,2 we performed extended haemodialysis (HD) treatment with high cut-(Theralite®, off (HCO) dialyzers Gambro) and bortezomib 1.3mg/m² + dexamethasone 20mg/day (days 1-4-8-11) therapy with improvement of renal function (creatinine 1,6mg/ dL) and normalization of FLC chains levels (k- and λ -FLC 5mg/dL and 6mg/ dL respectively with k/λ ratio=0.85) after the first chemotherapy cycle and 7 HCO-HD. Later, we performed 3 additional cycles with bortezomib dexamethasone at weekly administration (days 1-8-15-22 for a 35 days cycle) and a peripheral blood stem cells harvest for autologous transplantation program.

Conflict of interest

The authors declare that they have no conflicts of interest related to the contents of this article.

- Hutchison CA, Heyne N, Airia P, Schindler R, Zickler D, Cook M, et Al. Immunoglobulin free light chain levels and recovery from myeloma kidney on treatment with chemotherapy and high cut-off haemodialysis. Nephrol Dial Transplant 2012;27:3823-8.
- Hutchison CA, Cockwell P, Stringer S, Bradwell A, Cook M, Gertz MA, et al. Early reduction of serum-free light chains associates with renal recovery in myeloma kidney. J Am Soc Nephrol 2011;22:1129-36.

Gioacchino Li-Cavoli¹, Francesco Di Bassiano², Calogera Tortorici¹, Luisa Bono¹, Angelo Ferrantelli¹, Carlo Giammarresi¹, Rita Passantino³, Ugo Rotolo¹

- ¹ Division of Nephrology and Dialysis. Civic and Di Cristina Hospital. Palermo (Italy).
- ² Division of Ematology. Civic and Di Cristina Hospital. Palermo (Italy).
- ³ Anatomical Pathology Services. Civic and Di Cristina Hospital. Palermo (Italy).

Correspondence: Gioacchino Li Cavoli Division of Nephrology and Dialysis. Civic and Di Cristina Hospital, Palermo, Italy. gioacchinolicavoli@libero.it

Efficacy of High Permeability Haemodialysis in Acute Renal Failure due to Vancomycin

Nefrologia 2014;34(5):688-90

doi:10.3265/Nefrologia.pre2014.Jun.12630

To the Editor,

It was classically accepted vancomycin excretion through low permeability membranes during conventional haemodialysis (HD) was negligible (less than 5%). With high permeability membranes this has changed, which has important effects when treating vancomycin overdose. We report a case of acute oliguric renal failure secondary to vancomycin overdose resolved by the use of standard HD with high-permeability membranes.

An 83-year-old male with high blood pressure, type 2 diabetes mellitus, who had suffered a myocardial infarction 10 years before and with cognitive impairment due to Alzheimer's disease, was admitted with aspirative pneumonia. He was treated with amlodipine, metformin, aspirin, repaglinide and risperidone. On arrival he was suffering from moderate prerenal kidney failure (plasma creatinine [Cr]: 1.7mg/dl, estimated glomerular filtration rate calculated by the MDRD-4 formula $40 \text{ml/min}/1.73 \text{m}^2$), after being corrected stabilized (Cr: 0.95mg/dl, eGFR 80ml/min/1.73m²), with elemental anodyne urine. He was treated with furosemide piperacillin and tazobactam. Due to a poor prognosis, intravenous vancomycin was added on the fifth day at 0.5g/12 hours (13mg/kg/ day). Between the tenth and twelfth day urine volume decreased (200-400ml/ day) and Cr 6.3mg/dl (eGFR: 9 ml/ min/1.73m²) and elemental urine were seen suggestive of acute tubular necrosis (ATN). Doppler Ultrasound: preserved kidney size with no hydronephrosis with preserved intrarenal vasculature. Immunology and viral serology were negative. In 48 hours, Cr rose to 8.4mg/

dl (eGFR: 6.5ml/min/1.73m²). General condition worsened, with mild fever disorientation, drowsiness, dysarthria, and slight fever. Vancomycin levels (last dose 48 hours before): $87.7\mu g/ml$. No eosinophilia, eosinophiluria, skin or joint symptoms. The general condition of the patient was not appropriate for renal biopsy and the patient was diagnosed with acute kidney injury network (AKIN) III due to acute tubular necrosis secondary to vancomycin overdose. Acute HD was performed with a high permeability polyamide membrane (poliflux 21 H, Gambro®) 2.1m² surface, 50 microns thickness and ultrafiltration coefficient (KUF) 85ml/ hx mmHg for 4 hours, with a blood flow of 350 ml/min and 800ml/min. Following the session, vancomycin levels went down to 69.9µg/ml, due to an elimination of 20.3% and t_u of 12.6 hours. Two hours later, levels remained stable (65.5µg/ml). After a second session 24 hours later using an on line technique levels were at 41.6mg/ ml (percentage of reduction: 40.48% and t½: 6.13 hours). Three sessions of standard HD were performed with reduction percentages between 30-56% (Figure 1). Significant rebound was

only seen between the fourth and fifth session, separated by 72 hours. The patient's general condition gradually improved, fever disappeared and he recovered neurologically. During the following days, urine output increased progressively.

DISCUSSION


Vancomycin has a half-life of 4-6 hours (54-180 hours in terminal renal failure¹), 55% protein binding and 1500 Da molecular weight. It is eliminated by glomerular filtration. There are two causes of acute renal failure due to vancomycin acute interstitial nephritis (AIN) and ATN due to direct toxicity proximal convoluted tubules (PCT)² due to impurities related to its synthesis ("Mississippi mud"). Most cases present aggravating factors3: concomitant treatment with aminoglycosides or diuretics, treatment during more than 21 days, high APACHE score (Acute Physiology and Chronic Health Evaluation), obesity, trough concentration> 28µg/ ml and doses >4g/day. Cuprophane and cellulose triacetate membranes do not eliminate substances with a molecular weight >500Da, so that the decrease in concentration using these membranes is considered negligible. High permeability membranes have increased removal to 50%^{4,5}. The most efficient techniques are thought to be charcoal haemoperfusion, haemodiafiltration and high permeability standard HD, without greater effectivity of any having been clearly demonstrated.

Using biopsy diagnosis, we found 11 cases of acute renal failure: 63% had acute interstitial nephritis (AIN) and 45% acute tubular nephritis (ATN)6-¹⁰, and 4 cases that did not undergo biopsy, required high permeability HD or charcoal haemoperfusion. Patient ages ranged from 8 to 79, vancomycin doses from 1-2g/day, 40-60mg/kg/day, during 2-7 days in cases of ATN and 2-30 days in cases of AIN, with trough levels between 45.8 to 420µg/l. 92.8% of ATN cases recovered the baseline GFR despite needing acute HD. In all cases, high permeability membranes were used that achieved reductions of 50-73%.

In our case, the patient required five sessions. In the first, purification efficiency was mild (20.23%); in the second, on line, it was significantly higher (40.48%), and during the others vancomycin concentration reduction ranged from 30-56%. Clinical, neurological and renal function evolution were satisfactory, recovery of volume of diuresis and GFR in two weeks. In view of this case, we conclude that standard HD with highly permeable membranes is an effective treatment for acute renal failure caused by secondary ATN due to vancomycin overdose.

Conflicts of interest

The authors declare that they have no conflicts of interest related to the contents of this article.






Figure 1. Evolution of vancomycin levels and diuresis. HD: haemodialysis.

Nefrologia 2014:34(5):675-92

letters to the editor

- degrees of renal function. Antimicrob Agents Chemother 1984;25:433-7.
- Shah-Khan F, Scheetz MH, Ghossein C. Biopsy-Proven Acute Tubular Necrosis due to Vancomycin Toxicity. Int J Nephrol 2011;2011:436856.
- Bunchman TE, Valentini RP, Gardner J, Mottes T, Kudelka T, Maxvold NJ. Treatment of vancomycin overdose using high-efficiency dialysis membranes. Pediatr Nephrol 1999;13(9):773-4.
- Foote EF, Dreitlein WB, Steward CA, Kapoian T, Walker JA, Sherman RA. Pharmacokinetics of vancomycin administered during high flux hemodialysis. Clin Nephrol 1998;50:51-5.
- Ulinski T, Deschênes G, Bensman A. Large-pore haemodialysis membranes: an efficient tool for rapid removal of vancomycin after accidental overdose. Nephrol Dial Transplant 2005;20(7): 1517-8.
- Wu CY, Wang JS, Chiou YH, Chen CY, Su YT. Biopsy proven acute tubular necrosis associated with vancomycin in a child: case report and literature review. Ren Fail 2007;29(8):1059-61.
- Wicklow BA, Ogborn MR, Gibson IW, Blydt-Hansen TD. Biopsy-proven acute tubular necrosis in a child attributed to vancomycin intoxication. Pediatr Nephrol 2006;21(8):1194-6.
- SokolH, Vigneau C, Maury E, Guidet B, Offenstadt G. Biopsy-proven anuric acute tubular necrosis associated with vancomycin and one dose of aminoside. Nephrol Dial Transplant 2004;19(7): 1921-2.
- Htike NL, Santoro J, Gilbert B, Elfenbein IB, Teehan G. Biopsy-proven vancomycinassociated interstitial nephritis and acute tubular necrosis. Clin Exp Nephrol 2012;16:320-4.
- Panzarino VM, Feldstein TJ, Kashtan CE. Charcoal hemoperfusion in a child with vancomycin overdose and chronic renal failure. Pediatr Nephrol 1998;12:63-4.

Juan A. Martín-Navarro, Vladimir Petkov-Stoyanov, María J. Gutiérrez-Sánchez

Sección de Nefrología. Hospital del Tajo. Aranjuez, Madrid.

Correspondence: Juan A. Martín Navarro Sección de Nefrología, Hospital del Tajo, Aranjuez, Madrid. (Spain). juanmartinnav@hotmail.com

An atypical case of anti-GBM antibody disease with renal function deterioration from normal to end stage renal disease

Nefrologia 2014;34(5):690-2

doi:10.3265/Nefrologia.pre2014.Jun.12580

Dear Editor,

Anti-glomerular basement membrane (GBM) antibody disease accounts for 20% of all rapidly progressive glomerulonephritis. Occasionally it has been diagnosed in patients with normal renal functions and those patients have favorable renal prognosis. Here,

we report a patient with anti-GBM antibody disease who presented with nephrotic range proteinuria, normal renal function and lack of pulmonary symptoms. Contrary to the previous data in literature; he developed end stage renal disease (ESRD) within three years despite appropriate treatment.

A 23-years-old non-smoker university student male patient presented with minimal bilateral edema in lower extremities that was present for two weeks. He had no gross hematuria, hemoptysis or other pulmonary symptoms. His past medical history was unremarkable. Laboratory evaluation revealed 7200mg/day proteinuria with normal renal functions (serum creatinine: 1mg/dL). Serum albumin




Figure 1. First renal biopsy from the patient.

Large glomeruli are notable with minimal tubulointerstitial alteration (panel A). Higher power shows hypercellular glomeruli with crescents (arrow) and fibrinoid necrosis (arrow-head) (panel B). Lineer IgG staining along glomerular basement membranes on immunofluorescence (panel C) (A: H&E x 40; B: PAMS x 200; C: Immunofluorescence, FITC conjugated anti-IgG Ab x 200).